FLINN SCIENTIFIC, INC. Material Safety Data Sheet (MSDS)

MSDS #: 553.00

Revision Date: September 27, 2010

SECTION 1 — CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Nitric Acid Solution (<0.5 M)

Flinn Scientific, Inc. P.O. Box 219 Batavia, IL 60510 (800) 452-1261

CHEMTREC Emergency Phone Number: (800) 424-9300

SECTION 2 — COMPOSITION, INFORMATION ON INGREDIENTS

Nitric Acid (7697-37-2) < 6.3%, and Water (7732-18-5) 93-99%.

Synonym: 1 Molar or less CAS#: None Established

SECTION 3 — HAZARDS IDENTIFICATION

Colorless, odorless, liquid.

Irritating to body tissues. Avoid all body tissue contact.

FLINN AT-A-GLANCE

Health-0 Flammability-0 Reactivity-1 Exposure-1 Storage-1

0 is low hazard, 3 is high hazard

SECTION 4 — FIRST AID MEASURES

Call a physician, seek medical attention for further treatment, observation and support after first aid.

Inhalation: Remove to fresh air at once. If breathing has stopped give artificial respiration immediately.

Eye: Immediately flush with fresh water for at least 15 minutes.

External: Wash continuously with fresh water and mild soap for at least 15 minutes.

Internal: Rinse out mouth and give large quantities of water. May also give gastric antacids such as Milk of Magnesia. Do not induce vomiting, Call a physician or poison control at once.

SECTION 5 — FIRE FIGHTING MEASURES

Nonflammable, noncombustible liquid.

NFPA CODE

When heated to decomposition, emits highly emits toxic fumes of NO_x and hydrogen nitrate.

None established

Fire Fighting Instructions: Use triclass, dry chemical fire extinguisher. Firefighters should wear PPE and

SCBA with full facepiece operated in positive pressure mode.

SECTION 6 — ACCIDENTAL RELEASE MEASURES

Restrict unprotected personnel from area. Ventilate area and contain spill with sand or other inorganic absorbent material, neutralize with sodium bicarbonate or calcium hydroxide and deposit in sealed bag or container. See Sections 8 and 13 for further information.

SECTION 7 — HANDLING AND STORAGE

Flinn Suggested Chemical Storage Pattern: Inorganic #3. Store with amides, nitrates and nitrites.

Store in a dedicated acid cabinet and away from any source of water; if an acid cabinet is not available, store in Flinn Saf-Cube. Never store with acetic acid.

SECTION 8 — EXPOSURE CONTROLS, PERSONAL PROTECTION

Avoid contact with eyes, skin, & clothing. Wear chemical splash goggles, chemical-resistant gloves, & chemical-resistant apron. Use ventilation to keep airborne concentrations below exposure limits. Always wear a NIOSH-approved respirator with proper cartridges or a positive pressure, air-supplied respirator when handling this material in emergency situations (spill or fire). Exposure guidelines: TWA 2 ppm, STEL 4 ppm (OSHA, NIOSH)

FLINN SCIENTIFIC, INC.

Material Safety Data Sheet

Nitric Acid Solution (<0.5 M)

MSDS #: 553.00

Revision Date: September 27, 2010

SECTION 9 — PHYSICAL AND CHEMICAL PROPERTIES

Clear, colorless, odorless liquid.

Solubility: Soluble in water.

Formula: HNO₃ Formula Weight: 63.01 pH 1.0 (0.1 M solution)

SECTION 10 — STABILITY AND REACTIVITY

Avoid contact with bases, alcohols, alkaline earth metals, galvanized iron, aluminum, organic materials, amines, and acetic acid. Corrodes steel and other metals.

Shelf life: Good.

SECTION 11 — TOXICOLOGICAL INFORMATION

Acute effects: Skin and eye irritant

Chronic effects: N.A. Target organs: N.A.

ORL-HUMAN LDL0 (Concentrated): 430 mg/kg IHL-RAT LC₅₀: (Concentrated): 244 ppm/0.5 hr

SKN-RBT LD₅₀: N.A.

N.A. = Not available, not all health aspects of this substance have been fully investigated.

SECTION 12 — ECOLOGICAL INFORMATION

Data not yet available.

SECTION 13 — DISPOSAL CONSIDERATIONS

Please consult with state and local regulations.

Flinn Suggested Disposal Method #24b is one option.

SECTION 14 — TRANSPORT INFORMATION

Shipping Name: Nitric acid Hazard Class: 8, Corrosive UN Number: UN2031 N/A = Not applicable

SECTION 15 — REGULATORY INFORMATION

Not listed.

SECTION 16 — OTHER INFORMATION

This Material Safety Data Sheet (MSDS) is for guidance and is based upon information and tests believed to be reliable. Flinn Scientific, Inc. makes no guarantee of the accuracy or completeness of the data and shall not be liable for any damages relating thereto. The data is offered solely for your consideration, investigation, and verification. The data should not be confused with local, state, federal or insurance mandates, regulations, or requirements and CONSTITUTE NO WARRANTY. Any use of this data and information must be determined by the science instructor to be in accordance with applicable local, state or federal laws and regulations. The conditions or methods of handling, storage, use and disposal of the product(s) described are beyond the control of Flinn Scientific, Inc. and may be beyond our knowledge. FOR THIS AND OTHER REASONS, WE DO NOT ASSUME RESPONSIBILITY AND EXPRESSLY DISCLAIM LIABILITY FOR LOSS, DAMAGE OR EXPENSE ARISING OUT OF OR IN ANY WAY CONNECTED WITH THE HANDLING, STORAGE, USE OR DISPOSAL OF THIS PRODUCT(S).

Consult your copy of the Flinn Science Catalog/Reference Manual for additional information about laboratory chemicals.